



PCRBIO SYSTEMS
simplifying research

www.pcrbio.com

Clara® Probe 1-Step Mix Hi-ROX

Product description

Clara® Probe 1-Step Mix offers reliable probe-based qPCR detection of both RNA and DNA target sequences. Provided in a one-tube format, this powerful RT-qPCR mix gives superior target amplification, in single or multiplex assays, even from highly dilute samples.

Clara® Probe 1-Step Mix is a 4x qPCR mix containing hot start Taq polymerase, dNTPs, MgCl₂, an enhanced version of UltraScript® Reverse Transcriptase, and our RiboShield® RNase inhibitor, providing a complete 1-step RT-qPCR mix. It is developed to work well with the full range of probe types, including TaqMan, Scorpions and molecular beacons.

Quality control

PCR Biosystems operates under an ISO 13485 certified Quality Management System. Our products are extensively tested and undergo a comprehensive, multi-step quality control process according to ISO 13485 standards, to ensure optimum performance, consistency and traceability.

| Pack size | 4x Clara® Probe 1-Step Mix Hi-ROX |
|-----------------|-----------------------------------|
| 200 reactions | 1 x 1 mL |
| 600 reactions | 3 x 1 mL |
| 1000 reactions | 5 x 1 mL |
| 10000 reactions | 1 x 50 mL |

Shipping and storage

On arrival the kit should be stored between -30 °C and -20 °C. Avoid prolonged exposure to light. If stored correctly, the kit will retain full activity until the indicated expiry date. Avoid exposure of the stock solution to frequent temperature changes and limit handling at room temperature to the necessary minimum.

Technical support

Scan or click the QR code for troubleshooting help and answers to frequently asked technical questions. If you require further technical support, please email technical@pcrbio.com with the following information:

- Amplicon size
- Reaction setup
- Cycling conditions
- Screen grabs of amplification traces and melting profile



TROUBLESHOOT



FAQS

Product Use: Unless we agree otherwise in writing, the Goods we supply are provided:

1. For research purposes only and you should not use or rely on the Goods for diagnostic purposes. If you wish to use the Goods in a regulatory approved medical device, please contact us so that we may consider this and discuss it further with you.
2. Subject to our standard terms and conditions found at <https://pcrbio.com/terms-conditions/>.

Important considerations

Instrument compatibility: Different qPCR instruments may require different levels of ROX passive reference for normalisation. Use our qPCR Selection Tool to determine which ROX concentration your instrument requires (<https://pcrbio.com/resources/qpcr-selection-tool/>).

Template: The mix can be used with RNA or DNA extracted by most commercial kits or standard extraction methods, provided the amount and quality of template are within an acceptable range. Addition of 2 to 5 μL volumes of sample will improve assay precision.

Primer design: For efficient amplification we recommend amplicon lengths between 80-200 bp and not exceeding 400 bp. Shorter amplicons allow for faster cycling. Primers should have an approximate T_m of ~ 60 °C using default Primer 3 settings (<https://bioinfo.ut.ee/primer3/>). To verify the best annealing temperature for your primers in our products, please visit <https://pcrbio.com/resources/tm-calculator/>. For TaqMan probes choose a probe close to the 5' primer, avoiding terminal guanosine residues.

Reaction setup

1. Before starting, thaw and briefly vortex the 4x Clara® Probe 1-Step Mix Hi-ROX
2. Prepare a master mix based on the following table:

| Reagent | 20 μL reaction | Final concentration | Notes |
|------------------------------------|-------------------------------------|--------------------------|---|
| 4x Clara® Probe 1-Step Mix | 5 μL | 1x | |
| Forward primer (10 μM) | 0.8-2 μL | 400 nM - 1 μM | See above for optimal primer design |
| Reverse primer (10 μM) | 0.8-2 μL | 400 nM - 1 μM | |
| Probe (10 μM) | 0.25-1 μL | 125 - 500 nM | |
| RNA or DNA Template | 2-5 μL | Variable | <100 ng cDNA, <1 μg genomic DNA, 1 μg -1 μg total RNA, >0.01 μg mRNA, 4 to 1×10^8 copies viral RNA |
| PCR grade dH_2O | Up to 20 μL final volume | | |

3. Program the instrument using the following conditions, acquiring data on the appropriate channel(s) for your chosen probe(s):

| Cycles | Temperature | Time | Notes |
|----------------------|----------------------------------|--|---|
| 1 <i>Optional</i> | 52 °C | 5-10 minutes singleplex 10-20 minutes multiplex | Reverse transcription. Required only for RNA templates. |
| 1 | 95 °C | 3 minutes | Polymerase activation and RTase inactivation |
| 40-50 | 95 °C 55 °C-65 °C | 5-15 seconds 20-30 seconds | Denaturation Anneal/Extension |
| Melt analysis | Refer to instrument instructions | | Optional melt profile analysis, available for hybridisation probes only |