



**PCRBIO SYSTEMS**  
simplifying research

[www.pcrbio.com](http://www.pcrbio.com)

## UltraScript® 2.0 cDNA Synthesis Kit

### Product description

UltraScript® 2.0 cDNA Synthesis Kit is designed for fast, reliable and unbiased cDNA synthesis from a wide range of RNA sample types. The kit contains all the required components for cDNA synthesis, including a combination of anchored oligo(dT)<sub>18</sub> and random hexamers optimised to produce high quality cDNA for use in real-time PCR applications.

The kit utilises UltraScript® 2.0 Reverse Transcriptase (RTase), a robust and highly thermostable modified MMLV reverse transcriptase engineered for superior cDNA synthesis speed, yield and representation. The RTase is blended with an advanced RNase inhibitor to prevent degradation of RNA by contaminating RNase.

UltraScript® 2.0 RTase is not inhibited by ribosomal and transfer RNAs making total RNA an ideal substrate. The kit can be used with 20 µg to 3.5 µg total RNA or oligo(dT) purified mRNA.

The 5x buffer contains enhancers, anchored oligo(dT)<sub>18</sub>, random hexamers, dNTPs and MgCl<sub>2</sub>. The relative concentrations of random hexamers and anchored oligo(dT)<sub>18</sub> have been optimised for unbiased cDNA synthesis for use in subsequent real-time PCR experiments.

### Quality control

PCR Biosystems operates under an ISO 13485 certified Quality Management System. Our products are extensively tested and undergo a comprehensive, multi-step quality control process according to ISO 13485 standards, to ensure optimum performance, consistency and traceability.

Pack size	5x cDNA Synthesis Mix	UltraScript® 2.0 for cDNA Synthesis (with RNase Inhibitor)
25 reactions	1 x 100 µL	1 x 25 µL
100 reactions	4 x 100 µL	1 x 100 µL

### Shipping and storage

On arrival the kit should be stored between -30 °C and -20 °C. If stored correctly, the kit will retain full activity until the indicated expiry date. Avoid exposure of the stock solution to frequent temperature changes and limit handling at room temperature to the necessary minimum. Do not store the mix once it is combined with the RTase.

### Technical support

Scan or click the QR codes for helpful cDNA synthesis tips and answers to frequently asked technical questions. For further technical support, please email [technical@pcrbio.com](mailto:technical@pcrbio.com) with the following information:

- Amplicon size
- Reaction setup
- Cycling conditions
- Screen grabs of amplification traces and melting profile



TIPS & TRICKS



FAQS

**Product Use:** Unless we agree otherwise in writing, the Goods we supply are provided:

1. For research purposes only and you should not use or rely on the Goods for diagnostic purposes. If you wish to use the Goods in a regulatory approved medical device, please contact us so that we may consider this and discuss it further with you.
2. Subject to our standard terms and conditions found at <https://pcrbio.com/terms-conditions/>.

## Important considerations

**5x cDNA Synthesis Mix:** Contains anchored oligo(dT)<sub>18</sub>, random hexamers, 15 mM MgCl<sub>2</sub>, 5 mM dNTPs, enhancers and stabilizers. It is not recommended to add further enhancers or MgCl<sub>2</sub> to the reaction. The buffer composition has been optimised to generate high yield, non-biased cDNA for downstream applications.

**Template:** Use 20 pg to 3.5 µg total RNA or oligo(dT) purified mRNA for accurate quantification. Additional RNA is not recommended for quantification, as total reverse transcription is not guaranteed. As concentrations of target sequences will vary, users are encouraged to perform a template titration to find the optimal concentration for their application.

**Incubation temperature:** We recommend incubating with a temperature of 50-55 °C for 10-30 minutes for most applications. Where regions of interest contain high secondary structure (>65% GC), incubation temperatures of up to 70 °C may be used, but this will reduce the activity of the enzyme and may result in less total cDNA. The same temperature should be used when comparing samples.

**PCR and qPCR setup:** UltraScript® 2.0 contains a high concentration of RTase. Because excess RTase can inhibit DNA polymerase activity, we strongly recommend diluting the cDNA 10x-100x when it is to be used in PCR or qPCR reactions, even for low copy number gene expression analysis.

## Reaction setup

1. Allow 5x cDNA Synthesis Mix to thaw, briefly vortex.
2. Prepare a master mix based on the following table. Insert reagents in sequence listed:

Reagent	20 µL reaction	Final concentration	Notes
5x cDNA Synthesis Mix	4 µL	1x	
UltraScript® 2.0 for cDNA Synthesis (with RNase inhibitor)	1 µL		Add before total RNA as RNase inhibitor is blended with RTase
20 pg to 3.5 µg total RNA or polyA enriched mRNA	X µL		
PCR grade dH <sub>2</sub> O	Up to 20 µL final volume		

## No RT control setup (recommended for qPCR)

3. Prepare a master mix based on the following table. Insert reagents in sequence listed:

Reagent	20 µL reaction	Final concentration	Notes
5x cDNA Synthesis Mix	4 µL	1x	
20 pg to 3.5 µg total RNA or polyA enriched mRNA	X µL		Use equal amount as in step 2
PCR grade dH <sub>2</sub> O	Up to 20 µL final volume		

## Incubation and enzyme denaturation

4. Incubate at 50-55 °C for 10-30 minutes. Longer incubation times increase cDNA yields
5. Incubate at 95 °C for 10 minutes to denature RTase.
6. If using the resulting cDNA in qPCR or PCR please refer to “Important considerations” above.