



**PCRBIO SYSTEMS**  
simplifying research

[www.pcrbio.com](http://www.pcrbio.com)

## VeriFi® Library Amplification Mix

### Product description

VeriFi® Library Amplification Mix is a superior proofreading polymerase mix with AptaLock™ hot start technology. This 2x ready mix is designed for library amplification in NGS workflows and GC-rich PCR applications where the ability to amplify difficult targets without bias is required.

VeriFi® Library Amplification Mix contains the highly processive VeriFi® Hot Start Polymerase, developed for robust and versatile high fidelity PCR. The mix contains buffer, dNTPs, MgCl<sub>2</sub>, and enhancers and has been optimised to minimise GC-dependant bias during amplification. This unique mix composition enables improved library amplification, allowing for the acquisition of superior quality NGS datasets with a higher number of discrete reads than similar high fidelity mixes.

VeriFi® Library Amplification Mix can also be used to amplify difficult templates with very high or low GC content, when other proofreading enzymes fail, and can be used in multiplex PCR assays.

### Quality control

PCR Biosystems operates under an ISO 13485 certified Quality Management System. Our products are extensively tested and undergo a comprehensive, multi-step quality control process according to ISO 13485 standards, to ensure optimum performance, consistency and traceability.

Pack size	2x VeriFi® Library Amplification Mix
50 x 50 µL reactions	1 x 1.25 mL
250 x 50 µL reactions	5 x 1.25 mL

### Shipping and storage

On arrival the kit should be stored between -30 °C and -20 °C. If stored correctly, the kit will retain full activity until the indicated expiry date. The kit can be stored at 4 °C for 1 month.

### Technical support

Scan or click the QR code for our primer T<sub>m</sub> calculator and answers to frequently asked technical questions. For further technical support, please email [technical@pcrbio.com](mailto:technical@pcrbio.com) with the following information:

- Amplicon size
- Reaction setup
- Cycling conditions
- Screen grabs of amplification traces and melting profile



T<sub>m</sub> CALCULATOR



FAQS

**Product Use:** Unless we agree otherwise in writing, the Goods we supply are provided:

1. For research purposes only and you should not use or rely on the Goods for diagnostic purposes. If you wish to use the Goods in a regulatory approved medical device, please contact us so that we may consider this and discuss it further with you.
2. Subject to our standard terms and conditions found at <https://pcrbio.com/terms-conditions/>.

## Important considerations

**2x VeriFi® Library Amplification Mix:** The 2x mix contains VeriFi® Hot Start Polymerase, 6 mM MgCl<sub>2</sub>, 2 mM dNTPs, enhancers and stabilisers. Adding more PCR enhancers or MgCl<sub>2</sub> to the reaction is not recommended. The mix composition has been optimised to maximise PCR success rates.

**Primers:** For NGS library amplification, primers targeting the ligated adapters (e.g., P5 and P7 for Illumina platforms) should be used at a concentration between 0.4 µM and 1 µM. For standard end-point PCR, primers should have a predicted melting temperature between 60 °C and 70 °C, using default Primer 3 settings (<http://bioinfo.ut.ee/primer3/>) and the final concentration in the reaction should be between 0.2 µM and 0.6 µM.

**Denaturation:** Denaturation should be performed at 95 °C. However, if the presence of high GC regions results in low yields, increasing the melting temperature to 98-100 °C can improve the amount of product.

**Annealing:** We recommend performing a temperature gradient to experimentally determine the optimal annealing temperature. Alternatively, we recommend a 60 °C annealing temperature then increase in 2 °C increments if non-specific products are present.

**Extension:** Optimal extension is achieved at 72 °C. The optimal extension time is dependent on amplicon length and complexity of template. 30 seconds per kilobase (kb) is recommended for most applications. Two-step cycling protocols may also be used with combined annealing and extension at 68-72 °C.

**Cycle number:** Amplifying an NGS library excessively can lead to undesirable artifacts and cause amplification bias. The number of amplification cycles should be kept to the minimum necessary in order to obtain yields sufficient for downstream processes. Yields between 250 ng and 1000 ng are typically sufficient for most NGS applications. Depending on the DNA input, 5-15 cycles are generally enough. For non-NGS applications, 25-35 cycles are preferable.

## Reaction setup

1. Prepare a master mix based on the following table:

Reagent	25 µL reaction	50 µL reaction	Final concentration	Notes
2x VeriFi® Library Amplification Mix	12.5 µL	25 µL	1x	
Forward primer (10 µM)	1 µL	2 µL	400 nM*	*See above for optimal primer concentration
Reverse primer (10 µM)	1 µL	2 µL	400 nM*	
Template DNA	<100 ng genomic DNA <5 ng less complex DNA	<200 ng genomic DNA <10 ng less complex DNA	variable	
PCR grade dH <sub>2</sub> O	Up to 25 µL final volume	Up to 50 µL final volume		

2. Cycle using conditions based on the following table:

3-step cycling:

Cycles	Temperature	Time	Notes
1	95 °C	1 min	Initial denaturation
5-15 (for NGS), 25-35 (standard)	95 °C 55 - 68 °C 72 °C	15 seconds 15 seconds 30 seconds/kb	Denaturation (see "Important considerations" above) Anneal Extension

2-step cycling:

Cycles	Temperature	Time	Notes
1	95 °C	1 min	Initial denaturation
5-15 (for NGS) 25-35 (standard)	95 °C 68 - 72 °C	15 seconds 30 seconds/kb	Denaturation (see "Important considerations" above) Extension