

# qPCR BIO SyGreen Mix

- Sensitive
- Specific
- Fast

## Features

- Non-PCR inhibiting intercalating dye, better signal
- Rapid extension rate for early Ct values
- Market leading sensitivity - increased limit of detection
- Compatible on all real-time PCR platforms - standard and fast cycling conditions
- Blue mix available for easy sample visualisation during pipetting

## Applications

- Absolute quantification
- Relative gene expression analysis
- High throughput qPCR from genomic, cDNA and viral sequences
- Low copy number target genes

## Further Applications

- Crude sample PCR
- Standard and fast PCR conditions
- Specific amplification from complex templates (eg GC/AT rich)
- Compatible with all real-time PCR instruments

PCR Biosystems use a proprietary intercalating dye that does not inhibit PCR, unlike other popular fluorescent dyes. Combined with advanced enzyme, hot start and reaction buffer technology we offer market-leading sensitivity and reproducibility.

qPCR BIO SyGreen Mix can be used to quantify any DNA template including genomic, cDNA and viral sequences. Extremely low copy number targets can be detected specifically and with high efficiency. Antibody-mediated hot start technology prevents the formation of primer dimers and non-specific products leading to improved reaction sensitivity and specificity. Combining the latest developments in polymerase technology and advanced buffer chemistry we offer market-leading performance with minimal or no optimisation.




PCRBIOSYSTEMS  
simplifying research



# PCRBIOSYSTEMS

simplifying research



Black trace = Competitor Mix  
Purple trace = qPCR Bio SyGreen Mix

Figure 1.

Shows amplification and melt traces of 4 mouse housekeeping genes from a cDNA dilution series. qPCR Bio SyGreen Mix traces (purple) and 6 competitor mixes (black). Cycling conditions were 95°C 2min, 40 cycles of 95°C 10sec, 60°C 15sec on Roche LC480. For ACTG1 amplicon qPCR Bio mix was 2 to 4 Ct values earlier than 5 of 6 competitor mixes. The Ct was equal to that of Kapa Biosystems. The sensitivity of qPCR Bio mix was equal to 5 of 6 competitor mixes, but superior to Kapa Biosystems, demonstrated by absence of primer dimer at low template concentrations. For GAPDH amplicon qPCR Bio mix was 1 to 3 Ct values earlier for 4 of 6 competitor mixes and equal to 2 mixes. The sensitivity of qPCR Bio mix was superior to 4 of 5 competitor mixes, demonstrated by absence of primer dimer. Applied Biosystems mix showed equal sensitivity for this amplicon. For PGK amplicon, qPCR Bio mix had Ct values equal or lower than 5 of 6 competitor mixes. Sensitivity was equal to 4 mixes and superior to 2 mixes. For CANX amplicon, Ct values were 1 to 6 lower than 5 of 6 competitor mixes and equal to Kapa Biosystems mix. Sensitivity was superior to 3 of 6 mixes and equal to the other 3 mixes.

Overall, qPCR Bio SyGreen Mix outperformed each competitor mix on the 4 amplicons tested.

Catalogue Number	Product Name	Pack Size	Presentation
PB20.11-01	qPCR Bio SyGreen Mix Lo-ROX	100 x 20µl rxns	1 x 1ml
PB20.11-05		500 x 20µl rxns	5 x 1ml
PB20.11-20		2000 x 20µl rxns	20 x 1ml
PB20.12-01	qPCR Bio SyGreen Mix Hi-ROX	100 x 20µl rxns	1 x 1ml
PB20.12-05		500 x 20µl rxns	5 x 1ml
PB20.12-20		2000 x 20µl rxns	20 x 1ml
PB20.13-01	qPCR Bio SyGreen Mix with Fluorescein	100 x 20µl rxns	1 x 1ml
PB20.13-05		500 x 20µl rxns	5 x 1ml
PB20.13-20		2000 x 20µl rxns	20 x 1ml
PB20.14-01	qPCR Bio SyGreen Mix Separate-ROX	100 x 20µl rxns	[1 x 1ml mix] & [1 x 200µl ROX]
PB20.14-05		500 x 20µl rxns	[5 x 1ml mix] & [1 x 200µl ROX]
PB20.14-20		2000 x 20µl rxns	[20 x 1ml mix] & [4 x 200µl ROX]